



The Southeastern Biosafety Laboratory Alabama Birmingham (SEBLAB)

A top tier biocontainment research facility to advance your BSL-3 and ABSL-3 research program


Francisco Dominguez, Sawanan Saitornuang, Douglas Fox, Joseph Palmer, Kirk Vantrease,
Robert Aldredge, Adam Tehranchi, and Sixto Leal

Southeastern Biosafety Laboratory Alabama Birmingham (SEBLAB), Birmingham, AL 35294 Email: seblab@uab.edu



UAB The University of
Alabama at Birmingham

IMAGING (IN VITRO)



Instrument	Illumination, Objectives	Fluorophores	Resolution, Speed	Pros	Cons	Sweet spot(s)
A. Zeiss Lattice Lightsheet 7 Microscope	<ul style="list-style-type: none">LED (BF), 488 nm, 561 nm, 640 nm13.3x NA 0.4 Illumination44.9x NA 1.0 Detection	<ul style="list-style-type: none">BF, FL (490-740 nm)	<ul style="list-style-type: none">Voxel 145 nm x 145 nm x 145 nmVolume: 3 Volumes/s @ approx. 300 µm x 50 µm x 20 µmPlane: 400 frames/s @ approx. 300 µm x 20 µm	<ul style="list-style-type: none">Very low phototoxicityFast repeated imagingArivis analysis softwareCO₂ controller	<ul style="list-style-type: none">3 colors, no DAPILimited magnificationData sets can be very large	<ul style="list-style-type: none">Dynamic sub-cellular processes that requirerepeated (fast or long-term) imagingCell monolayers, organoids up to 200 µm thickPhagocytosis, vesicle trafficking, immune synapse, sub-cellular localization
B. Agilent Cytation 10 Confocal plate reader with Biostack 4	<ul style="list-style-type: none">LED (BF), 405 nm, 470 nm, 520 nm, 640 nm1.25x, 4x, 10x, 20x, 40x, 60x (air, long working distance)	<ul style="list-style-type: none">BF, FL (DAPI, GFP, Cy5, Texas Red, RFP), LUX, UV-VISMany additional objectives and cube sets available	<ul style="list-style-type: none">40 µm and 60 µm spinning disks96 well plate:Imaging: 8 min 9 secondsSweep mode: 10 seconds	<ul style="list-style-type: none">Up to 4 colorsConfocal resolution, magnificationPlate loader for fixed cell applicationsGen 5 analysis softwareCO₂/O₂ controller	<ul style="list-style-type: none">Incomplete imaging of well edgesRelatively slow imaging of live cellsData sets can be very large	<ul style="list-style-type: none">Plate reader assays (live or fixed)Long-term imaging with confocal resolution
C. Echo Revolve Fluorescence Microscope	<ul style="list-style-type: none">LED (BF), 405 nm, 488 nm, 555 nm, 640 nm4x, 10x, 20x, 40x, 100x (oil)	<ul style="list-style-type: none">BF, FL (DAPI, GFP, Cy5, Texas Red)	<ul style="list-style-type: none">3 channels, 111x25 µm square images in 5 seconds	<ul style="list-style-type: none">Simple, touch screen control, suitable for use with PAPR hoodFast imaging	<ul style="list-style-type: none">No long-term live imagingNo CO₂ controllerLimited analysis	<ul style="list-style-type: none">Multicolor, fast imaging without fixation


CORE RESEARCH TEAM

Available now to help PIs and research teams integrate new equipment into current and future research plans




PILOT GRANTS AVAILABLE


Starting 2026



MILTENYI MACSQUANT TYTO CELL SORTER AND FLOW CYTOMETER (BSL3)



IMAGING (IN VIVO)



Instrument	Location	Imaging Modes, Features, Animals	Pros	Cons	Sweet spot(s)
D. IVIS Lumina III	143 ABLS3	<ul style="list-style-type: none">2D BF, FL (red, far red), LUX<i>In vivo</i>, <i>ex vivo</i>Isoflurane anesthesiaMouse only, up to 5	<ul style="list-style-type: none">Repeated measurementsGood LUX sensitivity	<ul style="list-style-type: none">FL limited tissue penetrationFar red fluorophores are uncommonNot high resolution	<ul style="list-style-type: none">Track disease progression, dissemination, resolution
BIOEMTECH E. Beta Eye F. Gama Eye Imagers Scanners	143 ABLS3	<ul style="list-style-type: none">2D PET, BF2D SPECT, BF<i>In vivo</i>, <i>ex vivo</i>Isoflurane anesthesiaMouse only, 1 only	<ul style="list-style-type: none">Simple to useSensitive radioactive probesWide range of radionuclides can be used	<ul style="list-style-type: none">Combined pathogen, radioactivity hazardous materialsLimited commercial availability of probesOne animal at a time	<ul style="list-style-type: none">Track disease progression, dissemination, resolutionTrack drug, Ab probe localization, accumulation, elimination
G. MiLabs µCT	143 ABLS3	<ul style="list-style-type: none">3D µCTIsoflurane anesthesiaMouse, Rat, Hamster,	<ul style="list-style-type: none">Non-invasive, repeated imaging	<ul style="list-style-type: none">Limited scope of contrast agentsOne animal at a timeAnimals are exposed to x-rays	<ul style="list-style-type: none">Longitudinal tracking of clinically significant lesions

AEROBIOLOGY (ABSL3)



- Buxco Inhalation Exposure Tower
- Controlled nose-only exposure of pathogen, antigen or drug to mouse, rat, guinea pig, or ferret



- SciReq FlexiVent
- In vivo* lung function measurement for mouse, rat or ferret

COSMX SPATIAL MOLECULAR IMAGER



- Spatial single-cell transcriptomics
- FFPE or FF samples (BSL2)
- Human (HS), Mouse (MM)
- RNA Probe sets: 1k-plex (HS, MM), 6k-plex (HS), whole-transcriptome (HS), custom panels (HS, MM)
- RNA+protein multiomics (HS)
- Single-cell resolution (50 nm x,y)

CHROMIUM X CONTROLLER (AND BIORAD F100 THERMOCYCLER)



- GEM generation and cDNA preparation for 10X Genomics scRNAseq workflow
- Deployable to any BSL3 suite

A.



B.



C.



D.



Pandemic preparedness core (BSL2). A. ProMega Maxwell RSC 48 nucleic acid extraction system. B. Applied Biosystems QuantStudio 7 Dx Pro RT-PCR. C. Qiagen QIAcuity 4 digital PCR. D. Biotechne Protein Simple Ella Next Generation Elisa system.