

To make sure UAB is compliant with USDOT regulations, please be aware that only biohazard bags that meet ASTM-D Standards are appropriate for the disposal and transport of biohazardous waste.

Compliant Biohazard Bags-Primary Containment

Please see below for a list of biohazard bags that meet ASTM D1922 and impact resistance ASTM D1709 Standards.

Size	Manufacturer	Vendor	Quantity	Catalog Number
1 Gallon (8" x 12")	VWR	VWR	400/cs	14230-910
1 Gallon (8" x 12")	VWR	VWR	400/ca	14220-106
1.6 Gallon (11" x 14")	Medical Action	Amazon	50/box	50-42
1.6 Gallon (11"x14")	Medegen	Medline	500/cs	50-40
1.6 Gallon (11"x14")	Medegen	Medline	500/cs	50-42
3 Gallon (12" x 24")	VWR	VWR	200/cs	14220-090
5 Gallon (19" x 23")	Medegen/Fisherbrand	Fisher Scientific	200/cs	14828244
5 Gallon (19" x 23")	Medegen	Medline	400/cs	2634
5 Gallon (19" x 24")	VWR	VWR	200/cs	14230-912
5 Gallon (19" x 24")	VWR	VWR	200/cs	14220-110
10 Gallon (23" x 23")	Stericycle	Stericycle	200/cs	KP2323
10 Gallon (23" x 23")	Medegen	Medline	500/cs	F116
10 Gallon (23" x 23")	Medegen	Medline	400/cs	F116BX
10 Gallon (24" x 23")	Uline	Uline	200/cs	S-12984
10 Gallon (24" x 24")	Southland Medical	Southland Medical	250/cs	SMC122424
10 Gallon (24" x 24")	Medegen	Medline	150/cs	37-97
10 Gallon (24" x 24")	Medegen	Medline	500/cs	47-50
10 Gallon (24" x 30")	Medegen/Fisherbrand	Fisher Scientific	200/cs	14828245
10 Gallon (24" x 30")	VWR	VWR	200/cs	14220-094
10 Gallon (28" x 23")	Stericycle	Stericycle	500/cs	1523M
15 Gallon (24"x 30")	Southland Medical	Southland Medical	250/cs	SMC122433
15 Gallon (24" x 30")	Medegen	Medline	500/cs	F124
16 Gallon (24" x 32")	Medegen	Medline	250/cs	45-54
16 Gallon (25" x 35")	Medegen/Fisherbrand	Fisher Scientific	200/cs	14828246
16 Gallon (25" x 35")	Medegen	Medline	200/cs	2635
16 Gallon (28" x 31")	Medegen	Medline	250/cs	45-30
16 Gallon (24" x 36")	VWR	VWR	200/cs	14220-096
20 Gallon (30" x 36")	VWR	VWR	200/cs	11215-816
20 Gallon (30" x 36")	Uline	Uline	200/cs	S-20849
20 Gallon (30" x 36")	Medegen	Medline	250/Cs	F145
22 Gallon (30.5"x41")	Medegen	Medline	250/cs	F117

22 Gallon (31" x 41")	Medegen	Medline	100/cs	47-71
30 Gallon (31" x 38")	Medegen/Fisherbrand	Fisher Scientific	200/cs	14828247
30 Gallon (33" x 39")	ULine	ULine	100/cs	S-12985R
30 Gallon (39" x 33")	Stericycle	Stericycle	250/cs	2143
33 Gallon (31" x 43")	Southland Medical	Southland Medical	200/cs	SMC123143
33 Gallon (31" x 43")	Southland Medical	Southland Medical	100/cs	NON303143
35 Gallon (33" x 39")	Medegen	Medline	250/cs	F118
40 Gallon (37" x 48")	Medegen/Fisherbrand	Fisher Scientific	200/cs	14828248
40 Gallon (37" x 50")	Medegen	Medline	150/cs	F130
40 Gallon (37" x 50")	Medegen	Medline	150/cs	F110
40 Gallon (38" x 45")	Medegen	Medline	200/cs	F109
40 Gallon (38"x 45")	Medegen	Medline	100/cs	F169
40 Gallon (38" x 48")	VWR	VWR	200/cs	14220-098
40 Gallon (40" x 45")	Stericycle	Stericycle	100/cs	LIN042
40 Gallon (40" x 45")	Medegen	Medline	250/cs	F149
40 Gallon (40" x 46")	Medegen	Medline	100/cs	F135
40 Gallon (40" x 46")	ULine	ULine	100/cs	S-12986R
40 Gallon (45" x 40")	Stericycle	Stericycle	200/cs	2346
40 Gallon (40" x 46")	Southland Medical	Southland Medical	50/cs	SMC304046
44 Gallon (40" x 46")	Southland Medical	Southland Medical	100/cs	SMC123848
44 Gallon (40" x 46")	Southland Medical	Southland Medical	60/cs	SMC204046
55 Gallons (40"x 55")	Southland Medical	Southland Medical	50/cs	SMC304055
55 Gallons (40"x 55")	Medegen	Medline	150/cs	45-45
60 Gallons (55"x 60")	Medegen	Medline	50/cs	45-75
23"x8"x41"	Medegen	Medline	100/cs	47-41
23"x17"x46"	Medegen	Medline	100/cs	47-43
23"x8"x41"	Medegen	Medline	100/cs	47-45